## PRÁCTICA I.10

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## EXERCISE I.10

# ANALYSIS OF RACEMIC IBUPROFEN AND (S)-IBUPROFEN BY THE FORMATION OF DIASTEREOISOMERIC DERIVATIVES

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*E-mail: madalena@ff.up.pt*

### INTRODUCTION

Ibuprofen (4-isobutyl- $\alpha$ -methylphenylacetic acid) is a nonsteroidal anti-inflammatory agent widely used as an over-the-counter drug. Though only (2*S*)-ibuprofen is biologically active, commercial ibuprofen is currently marketed as a racemic mixture. Slow enzymatic racemization at the C-2 ensures that 57 to 71 % of the ingested drug is eventually converted into its active form [1,2]. Thus, while racemic ibuprofen requires almost 40 min to take effect, the same dosage of the (2*S*)-enantiomer requires less than a third of that time. For this reason, as well as to avoid the possible side effects caused by the use of racemates, pharmaceutical companies have become increasingly interested in producing enantiomerically pure drugs [2].

### AIM

The objective of this work is to distinguish (*R/S*)-ibuprofen from (2*S*)-ibuprofen.

### RESULTS AND DISCUSSION

Since the enantiomers have identical physical properties in an achiral environment, conventional methods of analysis are not able to distinguish them. Thus, racemic

ibuprofen is first converted into a mixture of diastereoisomeric amides by acid activation with 1,1'-carbonyldiimidazole and followed by reaction with *S*-(-)- $\alpha$ -methylbenzylamine (Fig. 1) [2]. The products resulting from this synthesis are then analyzed by thin-layer chromatography (TLC) (silica gel 60 F<sub>254</sub>; chloroform:petroleum ether 60–80°, 7:3). Racemic ibuprofen will generate (*S,S*) and (*R,S*) diastereoisomeric amides with the *R<sub>f</sub>* values of 0.44 and 0.37, respectively, while (*2S*)-ibuprofen will generate a (*S,S*) amide with the *R<sub>f</sub>* value of 0.37. Consequently, this simple method will allow us to distinguish (*2S*)-ibuprofen from the racemic mixture.

## MATERIAL AND METHODS

### Materials

Precoated silica gel 60 F<sub>254</sub> (Merck) TLC plates were used for TLC analysis. (*2R,2S*)-ibuprofen was purchased from Sigma Chemical Co. (*2S*)-ibuprofen and *S*-(-)- $\alpha$ -methylbenzylamine were obtained from Aldrich, and 1,1'-carbonyldiimidazole was obtained from Fluka.

### Methods [2]

The sample (50 mg) is dissolved in 0.9 mL of chloroform in a round-bottom flask. 0.3 mL of a solution of 1,1'-carbonyldiimidazole (0.14 g/mL in chloroform) is added dropwise, and the solution is stirred for 5 min at room temperature. 31  $\mu$ L of (*S*)-(-)- $\alpha$ -methylbenzylamine is then added by an automatic pipette to the reaction mixture. The solution is stirred for an hour at room temperature. The reaction is quenched by an addition of 3 mL of hexane and then washed with 2  $\times$  4 mL of a solution of sodium bicarbonate (5 %), 2  $\times$  4 mL of 1 M hydrochloric acid, 2  $\times$  4 mL of a solution of sodium bicarbonate (5 %), and finally with 2  $\times$  4 mL of water. After drying the organic layer with anhydrous sodium sulfate, the solution is evaporated to dryness to obtain the crude amide mixture, which is analyzed by TLC (silica gel 60F254; chloroform:petroleum ether 60–80°, 7:3).


Fig. 1 Synthesis of diastereoisomeric amides from (2*R*,2*S*)-ibuprofen and (2*S*)-ibuprofen.

## REFERENCES

- (a) M. S. Lennard, G. T. Tucker, H. F. Woods. In *Comprehensive Medicinal Chemistry*; C. Hansch, P. G. Sammes, J. B. Taylor (Eds.) *Principles of Pharmacokinetics and Metabolism*, J. B. Taylor (Ed.), Pergamon Press, England **5**, 196 (1990); (b) E. J. D. Lee et al. *Br. J. Clin. Pharmacol.* **19**, 669 (1985).
- S. E. Sen, K. S. Anliker. *J. Chem. Educ.* **73**, 569 (1996).

**Madalena Pinto**

madalena@ff.up.pt

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable

materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. HONORINA CIDADE (hcidade@ff.up.pt) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.