PRÁTICA I.9

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

E-mail: pradora@farmacia.ufmg.br

INTRODUCTION

Naproxen (6-methoxy- α -methyl-2-naphthaleneacetic acid) is an anti-inflammatory drug that inhibits cyclooxygenase, the enzyme involved in the biosynthesis of prostaglandins.

Naproxen is one of a class of enantiomerically pure drugs. Only the (*S*)-isomer is active and safe to use, because the (*R*)-isomer is reported to be a liver toxin. Therefore, the (*S*)-enantiomer is the eutomer.

Naproxen is used in the acid form and as its sodium salt. The eutomer of naproxen acid is dextrorotatory (specific optical rotation = $[\alpha]_D + 66$) and the eutomer of sodium salt is levorotatory ($[\alpha]_D - 11$). This offers a striking example of the lack of obvious relationship between structure and optical rotation, since both acid and salt have the same configuration about the chiral center.

**(S)-6-methoxy- α -methyl-2-naphthaleneacetic acid****sodium (S)-6-methoxy- α -methyl-2-naphthaleneacetate**

Specific optical rotation of chiral drugs is the parameter that indicates the enantiomeric purity. It is defined as optical rotation, expressed in degrees, observed (or calculated) with reference to the specific length of polarimeter tube of 1 dm and concentration of 1 g of solute in 1 mL of solution, measured under stated conditions (25 °C, D line of sodium). The specific optical rotation of a solid in solution is calculated by application of the following formula:

$$[\alpha]_D = (\alpha \times 100)/(l \times c)$$

where α is the observed angular rotation, l is the length of the cell in decimeters, c is the concentration in grams *per* milliliter.

The objective of the practical experiment described here is determining the specific optical rotation, $[\alpha]_D$, of naproxen and its sodium salt, both *S*-isomers.

EXPERIMENTAL

Basic reference: D. Walsh, S. Koontz. *J. Chem. Educ.* **74**, 585 (1997).

- Dissolve 150 mg, accurately weighed, of naproxen in chloroform to make 5 mL.
- Measure the zero reading of the polarimeter with chloroform.
- Transfer the solution of naproxen to the polarimeter tube. Measure the optical rotation of the solution.
- Dissolve 300 mg, accurately weighed, of naproxen sodium salt in methanol to make 5 mL.
- Measure the zero reading of the polarimeter with methanol.
- Transfer the solution of naproxen sodium salt to the polarimeter tube. Measure the optical rotation of the solution.

CALCULATIONS

Calculate the specific optical rotation of naproxen and its sodium salt and compare with the literature data.

OBSERVATION

After the experiment, remove the solvents of solutions and recover the naproxen and the naproxen salt.

Maria Auxiliadôra Fontes Prado

doraprad@dedalus.lcc.ufmg.b

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. Maria Auxiliadôra Fontes Prado who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.