PRÁCTICA IV.2

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EXERCISE IV.2

DETERMINATION OF THE CYTOTOXIC EFFECTS OF EXTRACTS, FRACTIONS, OR SUBSTANCES, BY MEANS OF THE MTT TEST

Clemencia de Castro de Pardo

Faculty of Medicine, Fundación Universitaria San Martín, Bogotá, Colombia

E-mail: cdepardo@sanmartin.edu.co

Purpose of the assay

To use the MTT redox assay in order to determine the potential cytotoxic effect of a given agent on cell lines tumor or primary cultures from the normal cells.

MTT redox assay

The assay is based on the metabolic bromide reduction from 3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazole (MTT), that is produced by the mitochondrial enzyme succinate-dehydrogenase, turned to a blue-colored compound named formazan. The mitochondrial functionality of the treated cells is then determined. This method has been extensively used to measure cell proliferation and survival capacities.

The remaining living cells are proportional to the amount of formazan produced.

The method was first developed by Mosmann in 1983, and later modified by Francois Denizot and Rita Lang in 1986.

Materials

- Living materials
 - Tumor cell line cultures
 - Primary normal cell culture
- Extracts, compounds, or substances to be tested
- Reagents and solutions
 - 0.2 % Dimethyl sulfoxide (DMSO)
 - Trypsin/EDTA
 - Fetal calf serum
 - Penicillin-streptomycin
 - MTT
 - Buffer phosphate (both Ca and Mg free)
 - Isopropanol
- Equipment
 - Laminar flow cabinet

Cell incubator (37 °C, 5 % CO₂)
 Centrifuge
 Micropipettes
 ELISA reader
 Inverted microscope

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Methodology


- Dissociate the cells using Trypsin/EDTA, resuspend them in supplemented culture media (antibiotics plus fetal calf serum).
- To be able to determine the required number of cells per well, it is necessary to perform the assay without the substance to be analyzed. The optical density (OD) must be 1 or more than 1.
- Considering the last fact, plant the cells into the 96-well microtitre plate; in most of the cell lines, 50.000 cells are used per well.
- The cells are incubated for 24 h at 37 °C, 5 % CO₂ to allow the attachment to the surface of the well.
- The sample (substance or extract to be analyzed) is prepared in DMSO 0.2 % (nontoxic concentration), in a way that the desired final concentration is reached adding 20 µL to the well.
- Add 20 µL of the sample solution to every well.
- Incubate 24 h at 37 °C, 5 % CO₂ to allow the substance to act on all the cells, including the case where it acts only in a specific phase of the cell cycle
- Discard the supernatant.
- Add to each well 100 µL of serum and phenol red free media.
- Add to each well 50 µL of MTT in a free Ca-Mg buffer.
- Incubate 4 h to allow formazan crystal formation.
- Discard the supernatant.
- Add 100 µL of isopropanol.
- Keep at room temperature until formazan crystals are dissolved.
- Read 450 nm OD.

Cell viability percentage is obtained as follows:

$$\% \text{ Viability} = \frac{\text{OD treated cells}}{\text{OD control cells}} \times 100$$

Each assay must be performed three times. Negative controls (solvent alone), positive controls (known cytotoxic substance), and untreated cells (OD 1 or more than 1) must be included.

Comparative values are calculated on a 50 % growth inhibition basis (GI₅₀), on the specific inhibitor agent treated cells.

Bibliography

Mosmann T. Rapid colorimetric assay for cellular growth and survival: Application to proliferation and cytotoxicity assay. *J. Immunol. Methods* 1983; 65: 55–63.

Denizot F., Lang R. Rapid colorimetric assay for cell growth and survival, Modifications to the tetrazolium dye procedure giving improved sensitivity and reliability. *J. Immunol. Methods* 1986; 89: 271–277.

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Clemencia de Castro de Pardo

cdeparado@sanmartin.edu.co

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories.

The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. Clemencia de Castro de Pardo (*cdeparado@sanmartin.edu.co*) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

The Colony formation techniques represent no risk for the students. The MTT technique requires the handling of this substance with gloves and within an extraction chamber because the substance is toxic.

If your exercise involves any "special" risks, please inform the editor.