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EXERCISE V.6
BENZODIAZEPINES – PART II
QUANTIFICATION OF THE ACTIVE SUBSTANCE IN COMMERCIAL DRUGS

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E-mail: norma@qo.fcen.uba.ar

To carry out the objective of the practice, you will have two situations:

1. If you know which is the active principle, then you will perform its isolation (see PART I), and the subsequent quantification.
2. If you have an unknown drug, previous to the quantification, you will perform the isolation, characterization, and identification of the active principle.

Quantitative Determination

From the benzodiazepines isolated from the commercial sample (see PART I), a verification of the active principle dose will be made using quantitative spectroscopic methods. Previously, a calibration curve must be made.

Appropriate Wavelength Determination

Prior to measurement, a determination of the appropriate wavelength must be made. In this particular wavelength, the sample must show an appreciable absorption intensity, to detect a low sample concentration.

On the other hand, it is recommended that the chosen wavelength shows few variations in the absorption intensity on the matter of small $\Delta\lambda$ in this specific case, the results will be reproducible and the different measurements of the curve will best fit.

Prepare a solution of the standard of the corresponding sample (approximately 10^{-4} M) and perform a qualitative UV spectrum. Pick an appropriate wavelength and perform the calibration curve measuring at this wavelength. The curve must be made using the diluted standard with methanol, at different concentrations; the molar concentration may be approximately 10^{-4} M or lower than this value, to observe the Beer and Lambert law.

$$A = \epsilon bc = \log I_0/I$$

ϵ = extinction molar coefficient.

b = cell length [cm]

c = solute concentration [g/l]

I_0 = incident light intensity

I = transmitted light intensity

Active Principle Quantification

Quantitative determination of the active principle in the sample will be made on the basis of the extrapolation of UV absorption, respective to the calibration

curve. Dissolve the extracted active principle in methanol (50 mL, using a graduated pipette), and use it to make dilutions 1:10 and 1:100. Determine the concentration of the original sample with these measurements and calculate the milligrams of benzodiazepine present in the original tablet.

If none of the measurements are in the useful linear zone of the curve, you must perform new dilutions from the original solution.

Miriam A. Martins Alho, Mirta L. Fascio, and Norma B. D'Accorso

norma@qo.fcen.uba.ar

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. Dra. Norma B. D'Accorso (norma@qo.fcen.uba.ar) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.