## PRÁCTICA VI.5

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\*E-mail: olock@pucp.edu.pe


## 1. INTRODUCTION

*Ruta* genus (Rutaceae family) includes around 60 species; the name *Ruta* derives from the Greek word *ruomai*, which means "to restrain", related to its supposed aphrodisiac action.

Species of this genus are commonly known as "ruda" and, from these, *Ruta graveolens* L. is an old native medicinal plant from Southern Europe. Among its various applications, the infusion of leaves is used as digestive, nervous tranquilizer, for intestinal gases, for liver pains, as well as for eye and ear inflammations. Other uses found are: against "the air" (rub temple with the leaves), against the fever (extract the plant in hot alcoholic beverage like brandy and apply on the head). It is also used in "chamanismo" practices to throw away the "bad spirits" (put a handful on the door or the house), etc.

The "rudas" are characterized by their flavonoid content, mainly the rutin. Rutin shows up as pale yellow crystals that can gradually become dark on light exposure. It decomposes with effervescence at 214–215 °C, has a specific rotation  $[\alpha]_D^{23} + 13.82$  (ethanol), and reacts with a dilute ferric chloride solution, developing a green color.

Rutin has an important use in capillary fragility treatment, for this reason, it is being included in pharmaceutical formulations.

In this experiment, we describe the detection and separation of the rutin from the "ruda" by usual chromatography techniques such as thin-layer (TLC), column (CC), and high-

performance liquid chromatography (HPLC). The experiment is performed in two sessions, developing points 2.2 and 2.3 in the first session and 2.4 in the second one.

## 2. EXPERIMENTAL PROCEDURE

### 2.1 REAGENTS

- Rutin, >98 % (as standard)
- RP-18 Solid-phase extraction cartridge, 500 mg
- Silica gel for column chromatography (0.063–0.2 mm)
- Silica gel for TLC (chromatoplates)
- Ferric chloride 1 % in ethanol
- Methanol, analytical grade
- Dichloromethane, analytical grade
- Ethanol, analytical grade
- Methanol, HPLC grade
- Water, HPLC grade
- Orthophosphoric acid, analytical grade

**CAUTION: Exposure to dichloromethane vapors should be avoided, so its manipulation should be done in the fume hood.**

### 2.2 QUALITATIVE ANALYSIS

#### 2.2.1 Extract preparation

- a. Macerate 10 g of dried and ground leaves with enough methanol in a stoppered flask during 48 h with continuous agitation.
- b. Filter the solution and evaporate it to dryness, in a rotary evaporator. Weigh the dry extract (this solution is labelled as E) and obtain the yield.

#### 2.2.2 Analysis of E solution by TLC

Rutin standard is used for comparison

- Adsorbent: Silica gel 60 F<sub>254</sub>
- Dimension of the plate: 10 x 4 cm
- Volume: 5  $\mu$ L (10 mg of E is dissolved in 1 mL of methanol)

- Mobile phase: CH<sub>2</sub>Cl<sub>2</sub>:MeOH (100:15)
- Detection: UV 366 nm, and 1 % FeCl<sub>3</sub> in ethanol

### 2.3 ISOLATION OF RUTIN BY COLUMN CHROMATOGRAPHY (CC)

#### 2.3.1 Sample preparation

- Dissolve 200 mg of E in approximately 5 mL of methanol.
- Add 1 g of silica gel (0.063–0.2 mm) and move it continuously until it is evaporated to dryness. This solid will be applied at the top of the column.

#### 2.3.2 Chromatographic conditions

- Adsorbent: Silica gel (0.063–0.2 mm, 70–230 mesh)
- Extract/adsorbent proportion (weight): 1/30
- Solvent composition for elution:

proportion		volume
CH <sub>2</sub> Cl <sub>2</sub>	MeOH	(mL)
100	15	30
100	20	30
100	30	30

#### 2.3.3 Column elution and analysis of fractions

- Apply the elute systems described in 2.3.2, collecting 10 mL volume fractions.
- Analyze all the fractions by TLC, with CH<sub>2</sub>Cl<sub>2</sub>/MeOH (100/15) as developing system.
- Combine the fractions that contain rutin and evaporate to dryness.
- According to the TLC results, purify by preparative TLC, with CH<sub>2</sub>Cl<sub>2</sub>/MeOH (100/15) as developing system.
- Recrystallize the solid and obtain the product yield.

#### 2.3.4 Spectral analysis

Obtain IR and UV spectra to analyze and/or compare with a standard or published data.

## 2.4 QUANTITATIVE ANALYSIS OF RUTIN BY HIGH-PERFORMANCE LIQUID CHROMATOGRAPHY (HPLC)

### 2.4.1 Calibration curve

Prepare a calibration curve with methanolic solutions containing 5.0, 15.0, 25.0, 35.0, and 50.0  $\mu\text{g}$  of rutin / mL.

### 2.4.2 Extract preparation

- a. Place 200 mg of the dried and ground plant material in a 20-mL Erlenmeyer flask.
- b. Add 7 mL of methanol and shake continuously for 20 min at room temperature.
- c. Filter the solution directly into a 25-mL volumetric flask.
- d. Repeat twice steps (b) and (c).
- e. Bring up to volume with methanol (this solution is labelled as X).

### 2.4.3 Solid-phase extraction (SPE)

- a. Condition the SPE column (RP-18 stationary phase) with 1 mL of methanol and then 1 mL of water. DO NOT DRY THE COLUMN.
- b. Apply 150  $\mu\text{L}$  of X solution into the column until the whole sample comes into the column.
- c. Elute at a rate of 0.5 mL/min (approximately) with 1.0 mL of MeOH:H<sub>2</sub>O (1:1). Receive the eluate in an appropriate container. This solution is for HPLC analysis (this solution is labelled as Y).

### 2.4.4 HPLC Analysis

- Column: RP - 18 stationary phase, 5  $\mu\text{m}$ , 125 x 4.0 mm, with guardcolumn RP – 18
  - Mobile phase:
    - A: Acetonitrile, with 0.05 % orthophosphoric acid
    - B: Water, with 0.05 % orthophosphoric acid.
  - Gradient program: 5 % up to 60 % A in B for 15 min
  - Flow: 1.0 mL/min
  - Monitoring wavelength: 254 nm
  - Retention time: 11–12 min

- Volume injection: 10  $\mu$ L

#### 2.4.5 Results

- Calculate rutin content, R, on the Y solution ( $\mu$ g of rutin / mL), using the calibration curve.
- Express rutin content on the original sample ("Ruda") by the formula:

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**Olga Lock**

[olock@pucp.edu.pe](mailto:olock@pucp.edu.pe)

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. OLGA LOCK SING ([olock@pucp.edu.pe](mailto:olock@pucp.edu.pe)) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.